

Venting a G-tube


1. Wash your hands for 15 seconds



2. Collect Supplies (60 ml syringe)



3. Pull the plunger out of the syringe and clamp the extension



4. Attach the syringe to G-tube extension



5. Unclamp the extension



6. Allow the air to come up through the tube.



7. You may hear air come up like a burp or fluid may fill the syringe, which is okay. Just be sure to let it all drain back in.




8. Clamp extension and take off the syringe




Giving Medications in a G-Tubes


1. Clean your hands and collect supplies (water, medicine, and extension)



2. Prime the extension with water 3-10 mls (depends on length of extension) unless already connected and primed



3. Clamp and attach the extension to the G-tube



4. Attach the medicine to the extension, unclamp, push in all the medicine, clamp, and then flush with water to push the medicine all the way through the G-tube



Anschutz Medical Campus 13123 East 16th Ave. Aurora, CO 80045 | 800-624-6553 | childrenscolorado.org

Author: Gastroenterology | Approved by Patient Education Committee | Valid through 2022

The information presented is intended for educational purposes only. It is not intended to take the place of your personal doctor's advice and is not intended to diagnose, treat, cure or prevent any disease. The information should not be used in place of a visit, call or consultation or advice of your doctor or other health care provider.

Children's Hospital Colorado complies with applicable Federal civil rights laws and does not discriminate on the basis of race, color, national origin, age, disability, or sex. • ATENCIÓN: si habla español, tiene a su disposición servicios gratuitos de asistencia lingüística. Llame al 1-720-777-1234. • CHÚ Ý: Nếu bạn nói Tiếng Việt, có các dịch vụ hỗ trợ ngôn ngữ miễn phí dành cho bạn. Gọi số 1-720-777-1234.